


Monitoring System

Large Screen Display

**No Coding Technology** 







**Before meal and after** 



## **SPECIFICATIONS**

Measurement Range	20~600 mg/dL (1.1~33.3 mmol/L)
Specimen Type	Capillary Whole Blood
Alternative Site Testing	YES (palm, forearm, upper arm, calf, and thigh)
Test Time	6 seconds
Sample Volume	0.7 μL
Sensitivity in HCT Range	20~60%
Operating Conditions	10~40°C (50~104°F) / 10~85% R.H.
	Altitude ≦ 3,402 m (11,161 ft) above
	the sea level
Storage/Delivery Conditions	4~40°C (39~104°F) / 10~85% R.H.
	(strip included)
	-20~50°C (-4~122°F) / ≦90% R.H.
	(strip excluded)
Memory Capacity	500 test results with date and time
Average Calculation	7, 14, 28, 60 and 90 days
Downloadable	YES (Micro USB port)
Battery Type	Two AAA 1.5V x2
Delete Memory	Yes
Meter Dimensions	100 mm (L) x 55 mm (W) x 18 mm (H)
Expiration for test strip	Storage: 30 months
	In-Use: 180 days

Manufacturer: OK Biotech Co., Ltd. www.okbiotech.com

Distributed Globally By:

Singapore's Emerging Biomedical Company www.jalmedical.com An assoclated company of OK Biotech Co., Ltd.

