



iCare

## iCARE RAPID IMMUNOGRAPHIC ASSAY MULTI-INFECTIOUS COMBO TEST (HIV 1&2 . HBSAG . HCV. TP)

### Multi-Infectious Disease Test (HIV/HBsAg/HCV/TP)

Rapid, Immunochromatographic Assay for Qualitative Detection of Infectious Disease (HIV/HBsAg/HCV/TP) in human whole blood, serum or plasma

#### Human Immunodeficiency Virus (HIV)

- Sensitivity: 100%
- Specificity: 100%
- Accurate results in 15 minutes

#### Herpatitis B (HBsAg)

- Sensitivity: 100%
- Specificity: 99%
- Accurate results in 15 minutes

#### Herpatitis C (HCV)

- Sensitivity: 100%
- Specificity: >99%
- Accurate results in 15 minutes

#### Syphilis (TP)

- Sensitivity: 100%
- Specificity: 99%
- Accurate results in 15 minutes

### Test Procedures



Apply sample to sample well




Add diluent to diluent well



Read results in 15 minutes

### Interpreting Results



Positive

Negative

Invalid

### Smart Options For Whatever You Want ...



Blood donation screening



Antenatal care system



Routine examination

