


iCARE DOA 1 Panel (Configurable) Rapid Screen Test

JAL jalinnovation.com/icare-doa-1-panel-configurable-rapid-screen-test/

27/02/2013

Category Name : Drug of Abuse (DOA) Rapid Screen Test Kit

- Rapid, qualitative, competitive binding immunoassay for the determination of the configured drug and its metabolites in urine at or above the cutoff level of 50ng/mL.
- Not intended to monitor drug levels, but only to screen urine for the presence of configured drug and its metabolites
- Conformed to USA, European and International drug of abuse test standards
- A variety of design specifications for your choice
- High specificity
- Simply to use
- Results in 5 minute
- Store in room temperature
- Easy to carry and suitable for Point-of-Care Testing (POCT)
- Tests for the following chemicals: AMP, BAR, BZO, BUP, COC, MDMA, KAN, THC, MTD, OPI, OCD, PCP, TCA



iCARE Drug of Abuse 1 Panel (Configurable)

Category Name	DOA Rapid Screen Test
Model Name	Refer to catalogue
Packing	50T / 100T per box
Specimen	Urine
Format	Panel

iCARE Drug of Abuse 1 Panel (Configurable)

Type of Test	Amphetamine (AMP), Barbiturates (BAR), Benzodiazepines (BZO) Buprenorphine (BUP), Cocaine (COC), Ecstasy (MDMA), Ketamine (KAN), Marijuana (THC), Methadone (MTD) Opiate (OPI), Oxycodone (OCD), Phencyclidine (PCP) or Tricyclic Antidepressant (TCA)
Specificity	100%
Sensitivity	90%
MOQ	2,000
Shelf Life	24 months from date of production
Private Labelling	Yes. MOQ: 5,000
Certification	CE, FDA, 510K for most similar products

[PreviousNext](#)

Made by [JAL Innovation \(S\) Pte Ltd](#)