

# iCARE DOA 10 Panels Rapid Screen Test (Configurable)

JAL [jalinnovation.com/icare-doa-10-panels-rapid-screen-test-configurable/](http://jalinnovation.com/icare-doa-10-panels-rapid-screen-test-configurable/)

28/02/2013

**Category Name :** Drug of Abuse (DOA) Rapid Screen Test Kit

- Rapid, qualitative, competitive binding immunoassay for the determination of the configured drug and its metabolites in urine at or above the cutoff level of 50ng/mL.
- Not intended to monitor drug levels, but only to screen urine for the presence of configured drug and its metabolites
- Conformed to USA, European and International drug of abuse test standards
- A variety of design specifications for your choice
- High specificity
- Simply to use
- Results in 5 minute
- Store in room temperature
- Easy to carry and suitable for Point-of-Care Testing (POCT)
- Tests for the following chemicals: AMP, BAR, BZO, BUP, COC, MDMA, KAN, THC, MTD, OPI, OCD, PCP, TCA



## iCARE Drug of Abuse 10 Panels (Configurable)

|               |                       |
|---------------|-----------------------|
| Category Name | DOA Rapid Screen Test |
| Model Name    | Refer to catalogue    |
| Packing       | 50T / 100T per box    |
| Specimen      | Urine                 |
| Format        | Panel                 |

## iCARE Drug of Abuse 10 Panels (Configurable)

---

|                   |   |
|-------------------|---|
| Type of Test      | Combination of any 10 of:<br>Amphetamine (AMP),<br>Barbiturates (BAR),<br>Benzodiazepines (BZO)<br>Buprenorphine (BUP),<br>Cocaine (COC),<br>Ecstasy (MDMA),<br>Ketamine (KAN),<br>Marijuana (THC),<br>Methadone (MTD)<br>Opiate (OPI),<br>Oxycodone (OCD),<br>Phencyclidine (PCP),<br>Tricyclic Antidepressant (TCA) |
| Specificity       | 100%  |
| Sensitivity       | 90%   |
| MOQ               | 2,000   |
| Shelf Life        | 24 months from date of production   |
| Private Labelling | Yes.<br>MOQ: 5,000  |
| Certification     | CE, FDA, 510K for most similar products   |

---

[PreviousNext](#)

Made by [JAL Innovation \(S\) Pte Ltd](#)