


# iCARE DOA 3 Panels (Configurable) Rapid Screen Test

JAL [jalinnovation.com/icare-doa-3-panels-configurable-rapid-screen-test/](http://jalinnovation.com/icare-doa-3-panels-configurable-rapid-screen-test/)

27/02/2013

**Category Name** :Drug of Abuse (DOA) Rapid Screen Test Kit



- Rapid, qualitative, competitive binding immunoassay for the determination of the configured drug and its metabolites in urine at or above the cutoff level of 50ng/mL.
- Not intended to monitor drug levels, but only to screen urine for the presence of configured drug and its metabolites
- Conformed to USA, European and International drug of abuse test standards
- A variety of design specifications for your choice
- High specificity
- Simply to use
- Results in 5 minute
- Store in room temperature
- Easy to carry and suitable for Point-of-Care Testing (POCT)
- Tests for the following chemicals: AMP, BAR, BZO, BUP, COC, MDMA, KAN, THC, MTD, OPI, OCD, PCP, TCA

## iCARE Drug of Abuse 3 Panels (Configurable)

Category Name	DOA Rapid Screen Test
Model Name	Refer to catalogue
Packing	50T / 100T per box

### iCARE Drug of Abuse 3 Panels (Configurable)

Specimen	Urine
Format	Panel
Type of Test	Combination of any 3 of: Amphetamine (AMP), Barbiturates (BAR), Benzodiazepines (BZO) Buprenorphine (BUP), Cocaine (COC), Ecstasy (MDMA), Ketamine (KAN), Marijuana (THC), Methadone (MTD) Opiate (OPI), Oxycodone (OCD), Phencyclidine (PCP), Tricyclic Antidepressant (TCA)
Specificity	100%
Sensitivity	90%
MOQ	2,000
Shelf Life	24 months from date of production
Private Labelling	Yes. MOQ: 5,000
Certification	CE, FDA, 510K for most similar products

[PreviousNext](#)

Made by [JAL Innovation \(S\) Pte Ltd](#)