



iCARE One-Step Drug Abuse Screening Test

① Instant ② Noticeable ③ Trustful ④ Effective ⑤ Convenient

Rapid, visual chromatographic assay for detection of abused drugs

- ▶ Sensitivity: >99.8%
- ▶ Specificity: >99%
- ▶ Unique 13 drugs test profile for urine specimen (AMP, BAR, BZO, COC, KET, MET, MTD, MDMA, OPI, PCP, THC, BUP, TCA)
- ▶ Unique 7 drugs test profile for saliva specimen (AMP, COC, THC, MET, BZO, MDMA)
- ▶ Special designs for NIDA/SAMSHA detection limits
- ▶ Multiple independent clinical studies and publications (CE, FDA, 510k)
- ▶ Can be customized for different abused drugs (multi 2-10 testing panel for urine; multi 2-7 testing panel for saliva)

