




Drug Testing Cup



A collection cup with testing device, delivers:

Instant results

Convenience

Accuracy

Tests are available in configurations of five to 16 of the following drugs:

- Amphetamine
- Benzodiazepines
- Marijuana
- Methamphetamine
- Methylenedioxymethamphetamine
- Morphine
- Phencyclidine
- Ketamine
- Oxycodone
- Barbiturates
- Cocaine
- Methadone
- Opiate
- Tricyclic Antidepressants
- Buprenorphine
- Propoxyphene



SPECIMEN COLLECTION


Open sealed pouch for sample collection. Sample level must be above the minimum fill line labeled on the cup.



Verify Specimen Temperature (Acceptable Specimen: 90-100°F / 32-38°C)

Remove the privacy label, and read results within **5** minutes.

RESULT INFORMATION



Want to become a distributor for TCup & Qcup?

Please contact the following in your region for more information about becoming one of our valued distributors, we are also offering **OEM** service.

▪ Instant ▪ Noticeable ▪ Trustful ▪ Effective ▪ Convenient

DOA



One Step Drugs-of-Abuse Test

Rapid, Visual Chromatographic Assay for Detection of Abused Drugs

- *Unique 13 drugs test profile for urine specimen (AMP, BAR, BZO, COC, KET, MET, MTD, MDMA, OPI, PCP, THC, BUP, TCA)*
- *Unique 7 drugs test profile for saliva specimen (AMP, COC, OPI, THC, MET, BZO, MDMA)*
- *Special designs for NIDA/SAMSHA detection limits*
- *Sensitivity: >99.8%*
- *Specificity: >99%*
- *Multiple independent clinical studies and publications (CE, FDA 510K)*



CE

190 Middle Road
#19-05 Fortune Centre
Singapore 188979
Tel: +65 8368 0800
Fax: +65 66185220
www.jalinnovation.com

JAL Innovation