

iCARE

RAPID IMMUNOGRAPHIC ASSAY INFECTIOUS DISEASE



One-Step Anti-Dengue (IgM & IgG) Triline Test

Rapid immunochromatographic assay for the qualitative detection of IgM & IgG type antibody to Dengue antigen detection in human whole blood, serum or plasma.

- Two Tests are available: serum or plasma test; whole blood, serum or plasma test.
- Sensitivity: 100%
- Specificity: 100%
- Accurate results in 15 minutes
- 24 months shelf life
- Easy to interpret color bands
- Customized




Test Procedures



Apply sample to sample well

15 minutes, read results

Result Interpretation



iCARE One-Step Anti-Dengue (IgM & IgG) Tri-Line Test