



# RAPID IMMUNOGRAPHIC ASSAY INFECTIOUS DISEASE



# Dengue


## ONE STEP Anti-Dengue (IgM & IgG) Tri-line Test

Rapid immunochromatographic assay for the qualitative detection of IgM & IgG type antibody to Dengue antigen detection in human whole blood, serum or plasma.

- Two tests are available: serum or plasma test; whole blood, serum or plasma test.
- Sensitivity: 100%
- Specificity: 100%
- Accurate results in 15 minutes.
- 24 months shelf life
- Easy to interpret color bands
- Customized



### Test Procedures:



Apply sample to sample well



15minutes. Read results

### Result Interpretation:



Catalogue Number	Description	Pack Size
ITP71001 TC40	ONE STEP DENGUE-IgM/IgG(serum or plasma)	40 Test Kit
ITP71001 DS50	ONE STEP DENGUE-IgM/IgG(serum or plasma)	50 Test Kit
ITP71002 TC40	ONE STEP DENGUE-IgM/IgG(whole blood,serum or plasma)	40 Test Kit
ITP71002 DS50	ONE STEP DENGUE-IgM/IgG(whole blood,serum or plasma)	50 Test Kit



190 Middle Road  
#19-05 Fortune Centre  
Singapore 188979  
Tel: +65 8368 0800  
Fax: +65 66185220  
www.jalinnovation.com