

iCARE Gonorrhoea Rapid Screen Test

JAL jalinnovation.com/icare-gonorrhoea-rapid-screen-test/

Category Name : Sexually Transmitted Disease (STD) Rapid Screen Test Kit

The iCARE One Step Gonorrhoea rapid home Screen test (swab) is a rapid home test kit for direct detection of Gonorrhoea.

This product is for in vitro diagnostic use only

Gonorrhoea is a common human sexually transmitted infection caused by the bacterium *Neisseria gonorrhoeae*.

With our iCARE Gonorrhoea Rapid Screen Test kit you can detect Gonorrhoea on your own in 100% privacy. Our test kit is easy to use and the result can be read within 15 minutes



Specifications

Category Name	STD Rapid Test Kit
Catalog No	WF1002
Packing	Hospital version: 50T Field version: 1T OTC Packing
Specimen	Whole blood / Serum / Plasma
Type of Test	Antigen
Format	Hospital version: 1 well cassette, pipette, diluents (4 bottles) Field version: Buffer A, Buffer B, extraction tube, swab, 1 well cassette, pipette
Sample	10ul serum
Result	Within 5 minutes
Specificity	98%
Sensitivity	98.8%

Specifications

Shelf Life 24 months

Storage Store at 4 °C ~ 30 °C in the sealed
pouch up to the expiration date.

Made by [JAL Innovation \(S\) Pte Ltd](#)