

iCARE HIV 1&2 Rapid Screen Test (Saliva)

JAL jalinnovation.com/icare-hiv-12-antibody-saliva-rapid-screen-test/

Category Name : Sexually Transmitted Disease (STD) Rapid Screen Test Kit

The iCARE HIV 1&2 Antibody Saliva Rapid Screen Test is a chromatographic immunoassay for the detection of antibodies to HIV 1&2 in human saliva. HIV 1&2 antibodies binding protein are precoated onto membrane as capture reagents in the Test region (T). During the testing, specimen is allowed to react with the colloidal gold reagents which have been labeled with HIV 1&2 specific antigens. Antibodies to HIV 1&2, if present, a pink colored band will develop on the membrane in proportion to the amount of HIV 1&2 antibodies present in the saliva specimen. Absence of this pink colored band in the Test region (T) suggests a negative result. To serve as a procedural control, a pink colored band in the Control region © will always appear regardless of the presence of absence of antibodies to HIV 1&2 in human saliva.



Made by [JAL Innovation \(S\) Pte Ltd](#)