



HIV 1&2 RAPID SCREEN TEST

Whole blood / Serum / Plasma



FOR FIELD & HEALTH CARE PROFESSIONAL USE

- Sensitivity >99%
- Specificity >99%
- Panel Confirmation 100%
- Detect HIV 1&2 and HIV 1 subtype "O"
- Medical advancement: Anti HIV Triline Test & HIV 1&2 Oral Rapid Test
- Less than 1 min hands on time
- Multiple independent parties evaluation including WHO

GMP

ISO 13485



Product of JAL Innovation Singapore

JAL Innovation

www.jalinnovation.com



HIV 1&2 RAPID SCREEN TEST

Whole blood / Serum / Plasma

ANTI - HIV TEST (WHOLE BLOOD / SERUM / PLASMA)

- Detect HIV Ab in whole blood / serum / plasma
- Detect HIV 1 & 2, and HIV 1 subtype "O" in whole blood / serum / plasma
- Detect antibodies of all subtypes (IgM, IgG, IgA)

ANTI - HIV TEST (URINE)

- Detect HIV Ab in urine
- Urine sample, require only 10 µl
- The safest method. No risk of infection

TEST PROCEDURES:

For Sandwich method (whole blood) & Captive method (urine / whole blood / serum / plasma)

STEP 1
Apply specimen into the well



STEP 2
Add 2 drops diluent into well




For Sandwich method (serum / plasma)


Apply specimen into the well



For One Step Anti-HIV Test



For One Step Anti-HIV Tri-line Test



Read result at 15 minutes

Refer to product inserts for each assay time

* Please refer to instruction insert for detailed testing procedures

GMP

ISO 13485



Product of JAL Innovation Singapore

www.jalinnovation.com