


iCARE

RAPID IMMUNOGRAPHIC ASSAY INFECTIOUS DISEASE



One-Step Malaria (P.f/P/v) Test

Rapid, Visual Chromatographic Assays for Detection of Malaria P.falciparum Specific Hifidine Rich Protein-2 (P.f HRP-2) and P.vivix Specific pLDH in Whole Blood


- Sensitivity: 96.9%
- Specificity: 99.7%
- Accurate results in 15 minutes
- 24 months shelf life
- Easy to interpret color bands
- Customized



Test Procedures



Result Interpretation



iCARE One-Step Malaria (P.f/P.v) Test