

■ Instant ■ Noticeable ■ Trustful ■ Effective ■ Convenient

MALARIA



One Step Malaria(P.f/P.v)Test

Rapid, Visual Chromatographic Assays for Detection of Malaria P.falciparum Specific Hifidine Rich Protein-2(P.f HRP-2)And P.vivix Specific pLDH in whole blood.

- Sensitivity: 96.9%
- Specificity: 99.7%
- Accurate results in 15 minutes
- 24 months shelf life
- Easy to interpret color bands
- Customized



CE

JAL Innovation

190 Middle Road
#19-05 Fortune Centre
Singapore 188979
Tel: +65 8368 0800
Fax: +65 66185220
www.jalinnovation.com


Rapid Malaria(P.f) Test

- Sensitivity:96.9%
- Specificity:99.7%
- Accurate results in 15 minutes
- Detection of circulating P. falciparum specific histidine rich protein-2 (P.f HRP-2) in whole blood
- 24 months shelf life
- Easy to interpret color bands
- Customized


Rapid Malaria(P.f / P.v) Test

- Sensitivity:96.9%
- Specificity:99.7%
- Accurate results in 15 minutes
- Detection of P. falciparum specific histidine rich protein-2 (P.f HRP-2) and P. vivax specific pLDH in whole blood
- 24 months shelf life
- Easy to interpret color bands
- Customized

Device Test Procedures:



Add the specimen to the sample well



Add the diluent to diluent well



Interpret the test results at 15 minutes

Result Interpretation:



Malaria P.f
Positive



Malaria P.f / P.v
Positive



Negative



Invalid

Catalog Number

Description

Pack Size

ITP11001 TC40	One Step Malaria (p.f) Test (Whole blood)	40Test Kit
ITP11003 TC40	One Step Malaria (p.f / p.v) Test (Whole blood)	40Test Kit

More Rapid Diagnostic Products Available

HBVs,HIV,HCV,Drugs of Abuse ,Syphilis,PSA,AFP,Pregnancy(HCG), Ovulation(LH),Cardiac Troponin I,TB,Malaria,H.pylori,Dengue,TOXO.

5 Seconds Blood Glucose Monitoring System **NEW**

Blood Grouping Test **NEW**