

▪ Instant ▪ Noticeable ▪ Trustful ▪ Effective ▪ Convenient

HCG



One Step Pregnancy Test

Rapid, Visual Chromatographic Assay for Detection of Early Pregnancy

- Two tests available: urine only and serum, plasma or urine test
- Detection limit: 25 mIU/ml for urine test
10 mIU/ml for serum/plasma/urine test
- Accuracy: (99% for both urine test and serum/plasma/urine test)
- Specificity: No cross reaction with high level LH, FSH, and TSH or any other Antibiotics
- Multiple independent clinical studies (FDA510K CE₀₁₂₃)



CE₀₁₂₃

190 Middle Road
#19-05 Fortune Centre
Singapore 188979
Tel: +65 8368 0800
Fax: +65 66185220
www.jalinnovation.com

JAL Innovation

One Step Ovulation Test

- High degrees of accuracy and reproducibility
- Control band used as a reference for semi-quantification of results
- Results in 3-10 minutes
- Accuracy:99%
- Multiple independent clinical studies

Test Procedures




Apply sample to sample well




Read results

Result Interpretation

For one step pregnancy test



For one step Ovulation test



More Rapid Diagnostic Products Available

HBVs,HIV,HCV,Drugs of Abuse ,Syphilis,PSA,AFP,Pregnancy(HCG), Ovulation(LH),Cardiac Troponin I,TB,Malaria,H.pylori,Dengue,TOXO.

5 Seconds Blood Glucose Monitoring System **NEW**
Blood Grouping Test **NEW**