

iCARE Vibro Cholerae (VCO 139) Rapid Screen Test

JAL jalinnovation.com/icare-vibro-cholerae-vco-139-rapid-screen-test/

26/02/2013

Category Name : Infectious Diseases Rapid Screen Test Kit

The *iCARE Vibrio Cholerae Rapid Screen* Test Kit is a chromatographic immunoassay for the qualitative detection of *V. Cholera* O139 in human feces, rectal swab specimens or water.

The test device has a test line “T” and control line “C” on the membrane. The control line is used for procedural control. Control line should always appear if the test procedure is performed properly and the test reagents of control are working. A pink T line will be visible in the result window if enough *V. Cholerae* O139 antigens are present in the sample. The special *V. Cholerae* O139 antibodies are used in the test device as both capture and detector materials. The detection limit for the VCO139 Test Kit is 10^4 cfu/ml. Sample equal to or greater than 10^4 cfu/ml will be tested positive.

**Generic
Packing**

Made by [JAL Innovation \(S\) Pte Ltd](#)